



## CF-Q160 I CF-Q160L KOLOSKOP VIDEOKOLOSKOP

SKU: AD042024

4.999,00 € 4.444,00 €

zzgl. [Versandkosten](#)

**Zustand:** Einsatzfähig – Sehr guter Zustand

**Lagerort:** Hamburg



---

**Categories:** [Flexible Endoskope & Zubehör](#)

### PRODUKTBESCHREIBUNG

- Channel Width: 3.7
- Working Length: L: 1680
- Field of View: 140
- Direction of View: Forward viewing Depth of Field: 3-100 Outer Diameter: 12.8 Outer
- Diameter Insertion Tube: L: 12.8
- Max Angulation Up: 180
- Max Angulation Down: 180
- Max Angulation Right: 160
- Max Angulation Left: 160 Key Benefits

### ZUSÄTZLICHE INFORMATIONEN

**Weight**

5 kg



## GALLERIE

