



**GE CZB 6.5 MHZ PART NR. 2152424 FOR LOGIQ 200**

**SKU: I102014-01**

CZB Curved Array Transducer 6.5 MHz. For small parts neonatal and biopsy applications.

Part Number 2152424

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



---

**Categories:** [Ultraschall / Doppler / Sonden /Zubehör](#)

## PRODUKTBESCHREIBUNG

CZB Curved Array Transducer 6.5 MHz. For small parts neonatal and biopsy applications.

Part Number 2152424

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg



[www.mtc-hamburg.com](http://www.mtc-hamburg.com)  
[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile.            0049 176 205 800 86  
Tel / Homezone   0049 40 18 20 79 03  
Tel / Büro         0049 40 3 61 61 - 256  
Fax.                0049 40 3 61 61 - 257

## ZUSÄTZLICHE INFORMATIONEN



## GALLERIE

