



OLYMPUS CLK 3 KALTLICHTQUELLE UND INSFFLATOR

SKU: AA092014

355,81 €

zzgl. [Versandkosten](#)

150 Watt Halogen

150 Watts Halogen Light Source

Has a dual lamp system allowing instant switch-over, in the event of lamp failure

The built in air pump allows insufflation and lens irrigation and makes it adaptable to all Olympus® endoscopes

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



Categories: [Flexible Endoskope & Zubehör,](#)
[Lichtquellen & Zubehör](#)

PRODUKTBESCHREIBUNG

150 Watt Halogen



150 Watts Halogen Light Source

Has a dual lamp system allowing instant switch-over, in the event of lamp failure

The built in air pump allows insufflation and lens irrigation and makes it adaptable to all Olympus® endoscopes

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN

