



## OLYMPUS GIF-K10 GASTROSCOPE, GASTROSKOP

**SKU:** AP092013

Field Of View : 100°

Insertion Tube : IT Diameter 11.4mm,

Angulation : Up 180°, Down 100°, Right 100°, Left 100°

Length : Working Length 103.5cm

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

0049 40 18 20 79 03

0049 176 205 800 86

fax.0049 40 65 86 33 64



---

**Categories:** [Flexible Endoskope & Zubehör](#)

## **PRODUKTBESCHREIBUNG**

Field Of View : 100°

Insertion Tube : IT Diameter 11.4mm,

Angulation : Up 180°, Down 100°, Right 100°, Left 100°

Length : Working Length 103.5cm

### **Compatibility**

[Olympus CLV-10 Light Source](#)



Olympus CLV-160 Light Source  
Olympus CLV-180 Light Source  
Olympus CLV-U20 Light Source  
Olympus CLV-U40 Light Source  
Olympus CLK-3 Light Source  
Olympus CLK-4 Light Source

Vollfunktionsfähig  
Im guten gebrauchten Zustand  
Artikel-Standort. Hamburg

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)  
[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)  
0049 40 18 20 79 03  
0049 176 205 800 86  
fax.0049 40 65 86 33 64

## **ZUSÄTZLICHE INFORMATIONEN**



## GALLERIE

