



[PENTAX ED-3430K DUODENOSCOPE / DUDONOSKOP](#)

SKU: AI092013 - AJ092013

2.499,00 €

zzgl. [Versandkosten](#)

Field Of View 120°

Direction Of View : Side

Depth Of Field : 4 - 100 mm

Insertion Tube : IT Diameter 11.3mm ,

Angulation : Up 120°, Down 90 °, Right 110 °, Left 90 °

Length: Working Length 125cm, Total Length 157cm

Biopsy Channel: Channel Size 3.8 mm

Compatibility : EPK-700, EPK-1000

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

www.mtc-hamburg.com


info@mtc-hamburg.com

Mobile 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



Categories: [Flexible Endoskope & Zubehör](#)

PRODUKTBESCHREIBUNG

Field Of View 120°



Direction Of View : Side

Depth Of Field : 4 - 100 mm

Insertion Tube : IT Diameter 11.3mm ,

Angulation : Up 120°, Down 90°, Right 110°, Left 90°

Length: Working Length 125cm, Total Length 157cm

Biopsy Channel: Channel Size 3.8 mm

Compatibility : EPK-700, EPK-1000

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile

0049 176 205 800 86

Tel / Homezone

0049 40 18 20 79 03

Tel / Büro

0049 40 3 61 61 - 256

Fax.

0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN

Weight

18 kg



GALLERIE

