



## **PENTAX EG-3430K GASTROSKOP / GASTROSCOPE**

**SKU:** AM092013

Field Of View : 140°

Direction Of View : Forward

Depth Of Field : 4-100mm

Insertion Tube : IT Diameter 11.4mm,

Angulation : Up 210°, Down 120°, Right 120°, Left 120°

Length : Working Length 105 cm, Total Length 137.7cm

Biopsy Channel : Channel Size 3.8mm

gepflegt

Inkl.Koffer

Vollfunktionsfähig

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

0049 40 18 20 79 03

0049 176 205 800 86

Fax.0049 40 65 86 33 64



---

**Categories:** [Flexible Endoskope & Zubehör](#)



## PRODUKTBESCHREIBUNG

Field Of View : 140°

Direction Of View : Forward

Depth Of Field : 4-100mm

Insertion Tube : IT Diameter 11.4mm,

Angulation : Up 210°, Down 120 °, Right 120 °, Left 120 °

Length : Working Length 105 cm, Total Length 137.7cm

Biopsy Channel : Channel Size 3.8mm

gepflegt

Inkl.Koffer

Vollfunktionsfähig

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

0049 40 18 20 79 03

0049 176 205 800 86

Fax.0049 40 65 86 33 64

## ZUSÄTZLICHE INFORMATIONEN



## GALLERIE

