



## **SHANDON CONSUL AUTOMATED ROBOTIC COVERSLIPPER**

**SKU:** ER042020-V

**2.378,81 €**

zzgl. [Versandkosten](#)

Consul is a robotic, fully automated glass coverslipping instrument capable of up to 400 slides per hour."

220v, 50/60hz,

Zustand. Gebraucht

Artike-Standort. Hamburg

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



---

**Categories:** [Labor](#), [Laborequipment](#)

## **ZUSÄTZLICHE INFORMATIONEN**



## GALLERIE

