



SIEMENS SOMATOM EMOTION 16

SKU: CT001



Categories: [Röntgen & Durchleuchtung](#), [Röntgen und zubehör](#), [Röntgen und Zubehör](#)

ZUSÄTZLICHE INFORMATIONEN



GALLERIE



```
System Serial Number : 37199  
Tube Serial Number   : 551301571  
Tube Type            : MCT_171  
Number of Scans      : 80365  
Kws Since Tube Inst. : 4.023e+006  
Tube Scan Seconds    : 338479  
System Scan Seconds  : 777896
```

