



**TOSHIBA PLG-805S / PLG 805 S 8 MHZ, 70
MM LINEAR PROBE**

SKU: BQ022016

799,00 €

zzgl. [Versandkosten](#)

Toshiba Linear 6.0 - 10.0 MHz PLG-805S Schilddrüse
Mamma Small Part

Toshiba

PLG-805S

Linear 6.0-10.0 MHz

Anwendung : Schilddrüse Mamma Small Part

Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



Categories: [Ultraschall](#) / [Doppler](#) / [Sonden](#)
[/Zubehör](#)

PRODUKTBESCHREIBUNG

Toshiba Linear 6.0 - 10.0 MHz PLG-805S Schilddrüse Mamma Small Part

Toshiba

PLG-805S



Linear 6.0-10.0 MHz

Anwendung : Schilddrüse Mamma Small Part

Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN

GALLERIE

