



WELCH ALLYN PIC40 E BIPHASIC AED TRAINING DEFIBRILLATOR

SKU: Q112016 // AD112016

888,00 €

zzgl. [Versandkosten](#)

Vollfunktionsfähig
Im guten gebrauchten Zustand
Artikel-Standort. Hamburg

www.mtc-hamburg.com
info@mtc-hamburg.com
Mobile 0049 176 205 800 86
Tel / Homezone 0049 40 18 20 79 03
Tel / Büro 0049 40 3 61 61 - 256
Fax. 0049 40 3 61 61 - 257



Categories: [Defibrillatoren](#)

PRODUKTBESCHREIBUNG

COLOR SCREEN

EXTERNAL PADDLES SET AND PEDIATRIC ADAPTER

Medicine calculations

NON INVASIVE BLOOD PRESURE WITH CUFF AND HOSE

CO2 (CAPNOGRAPH FEATURE) OPTION INCLUDED NO CABLE OR TUBING INCLUDED

OXYGEN SATURATION WITH ORIGINAL NELLCOR SENSOR AND WELCH ALLYN EXTENSION CABLE

ECG WITH 3 LEADS CABLE

PACEMAKER/ AED WITH ADAPTER

PRINTER WITH PAPER

DEFIBRILATION/EXTERNAL PACING VERIFICATION DEVICE

BATTERY AND POWER ADAPTER INCLUDED (100-220 volts)



Vollfunktionsfähig
Im guten gebrauchten Zustand
Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN

Weight 15 kg



GALLERIE

