



## WHALEDENT BIOSONIC UC100 ULTRASONIC CLEANER

**SKU:** BU112014

**399,00 €**

Model UC100

- 115 V - 1.1 A - 60 Hz

- Inside Chamber approx 10.5" X 6" X 2.75"

Vollfunktionsfähig

Artikel-Standort. Hamburg

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)


[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



---

**Categories:** [Dental und Zahnheilkunde](#)

## PRODUKTBESCHREIBUNG

Model UC100

- 115 V - 1.1 A - 60 Hz

- Inside Chamber approx 10.5" X 6" X 2.75"

Vollfunktionsfähig

Artikel-Standort. Hamburg



[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257

## ZUSÄTZLICHE INFORMATIONEN



## GALLERIE

